

## **PHYS 360/460**

### **Experiment 20 – NMR Spectroscopy**

#### ***Additional Notes***

##### **Saving Data**

- Acquire a single, non-updating scan to be clear about the waveform being saved (“Single” Button)
- If signal-to-noise ratio (SNR) is too low, use signal averaging ( $\text{SNR} \propto \sqrt{\# \text{ of scans}}$ )
- Set file quality to reduced (still 3000 – 5000 data points) to make the saved file size reasonable (~100 kB)
- When setting DC and TC before saving data, make sure:
  - the signal of interest is in view, while making it as large as possible
  - for multi-pulse experiments (IR, SE) both the signal of interest and the indicator for the first pulse (normally at the trigger point) are on screen to allow for accurate measurement of  $\tau$
  - include part of the baseline for later processing
- Make notes on the experimental setup (experiment type, sample type, rough idea of  $\tau$ ) for each file name, as these are numbered sequentially by default

##### **Data Pre-processing**

- The files will require baseline subtraction before they can be fit
- It is sometimes beneficial to perform filtering as well